

MULTI-BASE OVERDENTURE RESTORATION WORKFLOW




MULTI-BASE BALL ABUTMENT




MULTI-BASE LOCATOR ABUTMENT




Multi-base Ball abutments fixed in correct position in mouth.



Multi-base Locator abutments fixed in correct position in mouth.



Pick-up plastic impression copings attached on top parts.




Pick-up plastic impression copings attached on top parts.



Impression is transferred in stone model.



Impression is transferred in stone model.



Lab fabricates total denture with retainers.



Lab fabricates total denture with retainers.